



Reference: PTXBC047453

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACS2; FACL6; LACS 6; LACS2; LACS5; long-chain-fatty-acid--CoA ligase 6; fatty-acid-Coenzyme A ligase, long-chain 6; long fatty acyl-CoA synthetase 2; acyl-CoA synthetase long-chain family member 6		

For research use only.

GENE DESCRIPTION	acyl-CoA synthetase long-chain family member 6
GENE ID	23305
ORIGIN SPECIES	Human
ACCESSIONS	BC047453
VECTOR	pDONR223
NCBI SYMBOL	ACSL6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]