

DEFB1-defensin, beta 1 Gene

Reference: PTXBC047677

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BD1; DEFB-1; DEFB101; HBD1; beta-defensin 1; BD-1; defensin beta 1		

For research use only.

GENE DESCRIPTION	defensin, beta 1
GENE ID	1672
ORIGIN SPECIES	Human
ACCESSIONS	BC047677
VECTOR	pDONR223
NCBI SYMBOL	DEFB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagaacttccctctctgtgtttactctctgttactttgtctgagatggcctcaggtggttaactttctcacaggccttggccacagatctgatcattacaattgcgtcagcagtgaggaggcaatgtctctattctgcctgccgatctttacaaaattcaaggcacctgttacagagggaaggccaagtgtgcaagtga