

PSMA2-proteasome (prosome, macropain) subunit, alpha type, 2 Gene

Reference: PTXBC047697

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HC3; PMSA2; PSC2; proteasome subunit alpha type-2; macropain subunit C3; multicatalytic endopeptidase complex subunit C3; proteasome (prosome, macropain) subunit, alpha type, 2; proteasome component C3; proteasome subunit HC3; proteasome subunit alpha 2		

For research use only.

GENE DESCRIPTION	proteasome (prosome, macropain) subunit, alpha type, 2
GENE ID	5683
ORIGIN SPECIES	Human
ACCESSIONS	BC047697
VECTOR	pDONR223
NCBI SYMBOL	PSMA2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggagcgcgggtacagcttttgcgtgactacattcagcccgctctgtaaaactgtccagattgaatatgctttggcgtgctgtagctggaggagcccgctcgtgggaattaaagctgcaaatgggtgtggtattagcaactgagaaaaacagaaatccattctgtatgatgagcgaagtgt
acacaaagttagaaccaattaccaagcatatagggttggtgtacagtggtgcatgggccccgattacagagtgctgtgtcacagagctcgaaaactagctcaacaatactatctgtgtaccagaacccattctacagctcagctgtgtacagagagtagcttctgtgatgcaagaatatactcagt
caggtggtgtgtccattggagtttcttacttattgtgttggaatgaggagacacattatttattcagtcagatccatctggagcttactttgcctggaagctacagcaatgggaaagaactatgtgaatgggaagactttcttgagaaaagataatgaagatctggaacttgaagatg
ccattcatacagccatcttaaccctaaaggaaagcttgaaggcgaatgacagaggataacatagaagttggaatctgcaatgaagctggatttaggagcttactccaactgaagtaaggattacttgctgcatagcataa