

INO80E-INO80 complex subunit E Gene

Reference: PTXBC047712

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AI225782; AI854876; Ccdc85; Ccdc95; INO80 complex subunit E; NO80 complex subunit E; coiled-coil domain containing 95		

For research use only.

GENE DESCRIPTION	INO80 complex subunit E
GENE ID	283899
ORIGIN SPECIES	Human
ACCESSIONS	BC047712
VECTOR	pDONR223
NCBI SYMBOL	INO80E
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacggcggcgacggaagtggactacaaaaataaccggaatctgaagcggagctcaagttctcatctacgagcagcagtgcttcaggaggagctgaggaaagcgcaaggaaattactgaaggtgtccgggacaaagatttctcttagaccgacttctgcagt
acgagaacgtggatgaagactctcggaactcagatgccactgcatcatcagataacagcgagacggaggggacacccaagttgtctgacacacggcccttaagaggaaagagaagccctccgctggggggcgccccctctccctccagcctctccctgctcttcaacagggttccc
cttcaggcctccgggtccccccccatacctgagctcgtggtctctcccgctacccccattcctcttgactacctggccctgcagctgcccagagccagctccctgaggcccaagcgggagaaacggccccgctgccccgaaactcaagatggcgtgggacccccgact
gcccctgtggggggcgctgacctccctggccgggttctggggctggggctgggacacccctgacccccctccacccccctaagatgccccccccacgatcctgagcacggtccctcggcagatgttcagcgtgaggtagcgggacgatgcttggatggagacgatgacc
tgggtatcgacatcccgagtgga