

COX4I1-cytochrome c oxidase subunit IV isoform 1 Gene

Reference: PTXBC047869

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COX IV-1; COX4; COX4-1; COXIV; COXIV-1; cytochrome c oxidase subunit 4 isoform 1, mitochondrial; cytochrome c oxidase polypeptide IV; cytochrome c oxidase subunit IV; cytochrome c oxidase subunit 4I1		

For research use only.

GENE DESCRIPTION	cytochrome c oxidase subunit IV isoform 1
GENE ID	1327
ORIGIN SPECIES	Human
ACCESSIONS	BC047869
VECTOR	pDONR223
NCBI SYMBOL	COX4I1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttggtaccagggtatttagcctagttggcaagcagcaattccacctctgtgtgtgtacgagctcatgaaagtgttgaagagcgaagactttcgtccagcttatatggatcggtgaccaccccttgcggaggtggcccatgtcaagcacctgtctgccagccagaaggcac
tgaaggagaaggagaaggcctctgtgagcagccctccatggatgagaaagtcgagtggtggtattga