



Reference: PTXBC048287

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SEC63 homolog, protein translocation regulator; SEC63-like protein; SEC63-like; SEC63 protein translocation regulator; SEC63 homolog; translocation protein SEC63 homolog; DNAJC23; ERdj2; PCLD2; PRO2507; SEC63L		

For research use only.

GENE DESCRIPTION	SEC63 homolog (S. cerevisiae)
GENE ID	11231
ORIGIN SPECIES	Human
ACCESSIONS	BC048287
VECTOR	pDONR223
NCBI SYMBOL	SEC63
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]