



Reference: PTXBC049368

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TYRO3 protein tyrosine kinase; tyrosine-protein kinase receptor TYRO3; BYK; Dtk; Etk-2; RSE; Rek; Sky; Tif; tyrosine-protein kinase DTK; tyrosine-protein kinase RSE; tyrosine-protein kinase SKY; tyrosine-protein kinase TIF; tyrosine-protein kinase byk		

For research use only.

GENE DESCRIPTION	TYRO3 protein tyrosine kinase
GENE ID	7301
ORIGIN SPECIES	Human
ACCESSIONS	BC049368
VECTOR	pDONR223
NCBI SYMBOL	TYRO3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr