

RAB17-RAB17, member RAS oncogene family Gene

Reference: PTXBC050426

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB17, member RAS oncogene family; ras-related protein Rab-17		
For research use only.			

GENE DESCRIPTION	RAB17, member RAS oncogene family
GENE ID	64284
ORIGIN SPECIES	Human
ACCESSIONS	BC050426
VECTOR	pDONR223
NCBI SYMBOL	RAB17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcacaggcacacagaccgccagccagggctgccccagccagccccgtgtttcaagctgtttctctgggaagtggctccgtgggaagtccagcttggctcttcgggtacgtgaagaacgacttcaagagtatctgcctacgggtgggtgtgcttctcacaaggtgtg
gatgtgggtgccacctctgaagcttgagatctgggacacagctggccaggagaagtagcacagcgtctgccacctctactcaggggtgccaacgtgcgttctgtgtacgacatcaccaggaaggattcctcctcaaggctcagcagtggtgaaggacctggaggagctg
caccaggagaagtcctgtgtgctgtgtgggcaacaagacggacctcagccagagcgggaggtgacctccaggaagggaaggtttgccgacagccagaagttgctgttcatggaacttcggccaaactgaaccacaggtgtcggaggtgttcaatacagtggtccaagag
ctactgcagagaagcgacgaggaggccaggtctacgggggatgcagctgtgtgtgaacaaggggcccgaggcaggccaaatgtcgcgccactag