



Reference: PTXBC050454

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DRS2; DVL; DVL1L1; DVL1P1; segment polarity protein dishevelled homolog DVL-1; DSH homolog 1; dishevelled 1 (homologous to Drosophila dsh); dishevelled, dsh homolog 1; dishevelled-1; dishevelled segment polarity protein 1		

For research use only.

GENE DESCRIPTION	dishevelled, dsh homolog 1 (Drosophila)
GENE ID	1855
ORIGIN SPECIES	Human
ACCESSIONS	BC050454
VECTOR	pDONR223
NCBI SYMBOL	DVL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]