

LOC202051-hypothetical protein LOC202051 Gene

Reference: PTXBC050563

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC202051
GENE ID	202051
ORIGIN SPECIES	Human
ACCESSIONS	BC050563
VECTOR	pDONR223
NCBI SYMBOL	LOC202051
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgacgccctcgggtggtcgaaggcgggtcaggatctgtgtctcgtttagatcaactgcgggacgtgattgagtcctcaggaggaactaaccacagctgaggaaacgtgatgttctccaggacgaaattttgtagtaagaagagttccaggcagtggaagaagctg
gtggaagagaaagctgccatgccaaacaaaggtctcctgccaagggaaggagagagttacagtttccctcggagaggtagagggtctatccaagcagctggagaagagagaagctggcctttgaaaaagcgtctccagtgtaagagcaaatcttcaggagtcagcaag
aaggaccagctcatcaccaagtgcattggtaggcgggaggccgtgtctccaccacggttatataaacgggctctttcctaaccaccacagtgctgtgtcaaaaggcagggtttcacttaacacgtgtccacagctgaggatagcaggcccggtgggcccacaccttcgaaaggc
cttaaggaccaggcaggccacggaaatga