

YPEL3-yippee-like 3 (Drosophila) Gene

Reference: PTXBC050664

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein yippee-like 3; yippee like 3		
For research use only.			

GENE DESCRIPTION	yippee-like 3 (Drosophila)
GENE ID	83719
ORIGIN SPECIES	Human
ACCESSIONS	BC050664
VECTOR	pDONR223
NCBI SYMBOL	YPEL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaacctgaggcctctcactgagtcacaaacgagacgcctctggtecccgagagggctgctgctctgggaaaggagagagcagccagtggtgacagagcagtgccccacggctctgggagtgccctgcatgtgtgtggccaggctctgacagccacactactccctcccg
caggcactgggctccctctgctcccggtggccgctcccgctggggccactgccccggccccgccccccatgggtgctggatgcaagcccaagacgtttcaggcctactggatgattgtaccggagggtatagctgtgccactgcccgcctcactggccaaccacgacgacat
ctccaagtcttccagggcagtcagggcgctgacctcttcaactcagtggtgaacgtgggctgccccgagccgaggagcgggtgctgctgacggcctccatgctgctgacatccactgcgagaactgcaagaccatttgggctggaatatgaacaggccttggagagca
gccagaagtacaaagggggaagtacatcattgaactcaaccatgatcaagacaacggctgggactga