



Reference: PTXBC050666

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HMGA1L4; HMGIYL4; IES2; PAP-1BP; PAPA-1; PAPA1; ZNHIT4; hIes2; INO80 complex subunit B; IES2 homolog; PAP (Pim-1 associated protein) associated protein 1; PAP-1 binding protein; PAP-1-associated protein 1; high mobility group AT-hook 1-like 4; zinc finger		

For research use only.

GENE DESCRIPTION	INO80 complex subunit B
GENE ID	83444
ORIGIN SPECIES	Human
ACCESSIONS	BC050666
VECTOR	pDONR223
NCBI SYMBOL	INO80B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr