

RPSA-ribosomal protein SA Gene

Reference: PTXBC050688

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	37LRP; 67LR; ICAS; LAMBR; LAMR1; LBP; LBP/p40; LRP; LRP/LR; NEM/1CHD4; lamR; p40; 40S ribosomal protein SA; 37 kDa laminin receptor; 37/67 kDa laminin receptor; 67 kDa laminin receptor; colon carcinoma laminin-binding protein; laminin receptor 1 (67kD, ri		

For research use only.

GENE DESCRIPTION	ribosomal protein SA
GENE ID	3921
ORIGIN SPECIES	Human
ACCESSIONS	BC050688
VECTOR	pDONR223
NCBI SYMBOL	RPSA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccggagcccttgatgtcctgcaaatgaaggaggagatgtccttaagttccttcagcagggaacccacttaggtggccaactctgactccagatggaacagatcatctataaaaggaaaagtgtggcatctatatcataaatctcaaggagacctgggagaagcttctgtgcca
gctcgtgcaattgttgcattgaaaacccctgctgagtcagtgttatatctccaggaaactggccagaggcctgtgctgaagtgtgtgctgcccactggagccactccaaattgctggccgcttactctggaaccttactaaccagatccaggcaaccttccggagccacggcttctgt
ggttactgacccaggcgtgaccaccagcctctcacggagcattctatgttaacctactaccattgcgctgtgtaacacagattctcctctgctgtatgtggacattgccatccatgcaacaacaaggagctcactcagtggtttgatgtgggtgatgctggctcgggaagttctgcgc
atgcgtggcaccattccctgaacacccatgggaggtcagcctgactgtactctacagagatcctgaagagattgaaaaagaagagcaggctgctgctgagaaggcagtgaccaaggaggaattcagggtgaatggactgctccgctcctgagttcactgctactcagcctgaggt
tgcagactggtctgaaggtgtacaggtgccctctgtcctattcagcaattccctactgaagactggagcgtcagcctgcccaggaagactggctgctgagctccactgctcagggcactgaatgggtaggagcaaccactgactggtcttaa