

C4orf30-chromosome 4 open reading frame 30 Gene

Reference: PTXBC050697

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C4orf30; DDB1- and CUL4-associated factor 16; DDB1 and CUL4 associated factor 16		
For research use only.			

GENE DESCRIPTION	chromosome 4 open reading frame 30
GENE ID	54876
ORIGIN SPECIES	Human
ACCESSIONS	BC050697
VECTOR	pDONR223
NCBI SYMBOL	C4orf30
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtcctagaaatcctctcctgaccactgtcagaatcagaagtggaggaagaagaaatattagttacataaatgagagttctggggaagagtggaattcctctgaagaagaggactctatgggtcccaactatgcctcttgagagttcctggcagggttaagtgcctttaaatatt
ccacaacttggaacctttaaatcctaatcctggtgtatcatgctaaactgttgatccaagcacaccagtcctatacttcgagagataggctcgaagactctccattgtccattgtgtcccaactggaaccaattcctgaatggccctctggcctctgtggagtcacacctttcaaaa
gcctcttacaagtccagccgctctctagagatcatccacttaaatggagcactgcaatttgcaccaaacagcctaagccgaacattgagtagagccactcccatcctgaatacctaaacagatccctaattcatgtttcttgggtgtgtgtgctggtgactaaaacagtttaagg
aaacaactcgtactgaacctcaacactatttcttactgacttccaaaaggcagtttaacaaactcctaactgcacactgtaa