

KIAA0562-KIAA0562 Gene

Reference: PTXBC050721

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KIAA0562; CFAP256; GlyBP; JBTS25; ROC22; centrosomal protein of 104 kDa; centrosomal protein 104kDa; centrosomal protein 104		

For research use only.

GENE DESCRIPTION	KIAA0562
GENE ID	9731
ORIGIN SPECIES	Human
ACCESSIONS	BC050721
VECTOR	pDONR223
NCBI SYMBOL	KIAA0562
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccccacaagattggattttagtctcagctcatctggacacgaagacggcttcagtgccccggagctcatgatccacgcgccaactgtcagtggggtggcggtcacctagatttgcagttccacaagaattgtcctcaaatggaggagatgtcgaataaggaaactgcagttac
ttgctcaccagtatatgattcaagtaaaattgagttctacattagtgaaagcttgccctgaataatttgcaccctatcaagcagagcgggtttcgaagactggctacgtgtctctctgtgataatgaaaagacaggttgcaaacccgggaactaaatcagtttatgtggatgcagtaggacaatttc
ttaaactgattttcaccacaaacctgtcaacaatacaacatataatcaggttgcttgggtgccataataatcattggagacctgcagatttcagtgatgaaagcaatactgcctctcagagagaattgattgaccactaccttgggcacaacagcagggaccctgctctagaaggaaacgta
cgccagcagtgtaagaacccgggggtgaatccactttggggagctaaaggccccgctgtgcatcctcgtgcacctctctgtctgggcacctcttgggacagtggtccttgcacctgcctgctgtggatgacaatgaggggactccattccagaggtgcctgtctaa