

C20orf103-chromosome 20 open reading frame 103 Gene

Reference: PTXBC050727

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LAMP family protein C20orf103; C20orf103; BAD-LAMP; BADLAMP; LAMP-5; UNC-46; lysosome-associated membrane glycoprotein 5; brain and dendritic cell-associated LAMP; brain-associated LAMP-like protein; lysosome-associated membrane protein 5; lysosomal assoc		

For research use only.

GENE DESCRIPTION	chromosome 20 open reading frame 103
GENE ID	24141
ORIGIN SPECIES	Human
ACCESSIONS	BC050727
VECTOR	pDONR223
NCBI SYMBOL	C20orf103
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatctccaaggaagagggtcccccagcatcgacagactcgaattctctgatgttttccatacaatggtcacaatcatggcagaacaagaagtggaaaatctctcaggcccttcactaacctgaaaaagatatattgtgtgctgggaaaaatgggacgacgtgtctcatggcagag
tttgcagccaaatttattgtacattatgatgtgtggccagcaactacgtatgatcgcagacagccgatatcgcatggaccggggagctgaggtgaaggccgctgtggccacagccagtcggagctgcaagtgttctgggtgcatcgcatatgcactcaaatgctcttga
aagaaagccacaacatgtccaaggacgtgagcgacttgagcgctgagcaaatgcaattgttctacgactctcgggagaaaacccacttcaagacgcagtcagtgctgggaagcacacagccaactgcaccacctctctgcttgcaccccgctgggaagtctatgagtg
tcaagctcaacaacatttactggcctctagtatccgcagaagacggtcaccatgatcctgtctggtgccatccacacttttgacattatctcagattttgtctcagtgaaagacataaatgccagtgatgagcgggagcaactggaagaaccttgccttgattttgggctcat
cttggcctcgtcatcatggttaacactcgcatttaccacgtccaccacaaaatgactgccaaccagggtgcagatccctcgggacagatccagttataagcacatgggctag