

C11orf74-chromosome 11 open reading frame 74 Gene

Reference: PTXBC051031

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C11orf74; HEPIS; NWC; chromosome 11 open reading frame 74		

For research use only.

GENE DESCRIPTION	chromosome 11 open reading frame 74
GENE ID	119710
ORIGIN SPECIES	Human
ACCESSIONS	BC051031
VECTOR	pDONR223
NCBI SYMBOL	C11orf74
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgcccattatgcaggattggaataatggatgaagatcaatcaaaagacgtcttgataaattccttaattgtcatgagcaacatatgatgaagaattctgaacacttttactcatcttcacaaggagcatgtgtccaaagggaggtgttggaactgattctcagaaaacatttt
tacctcagcaaaagtactcataaaatgaagcagatgactaccattcttagaataaaacacattcttctgacttcacacaatgttggaagaacaggtagataatttctagattagaagatttggacatggatgaagagattaaacccaaatgagtgaggatttgccttcaggag
aagtggagcaggatgtaagcaccagcatctctctgtatcccttttggccagcctcctacgtgtgaagtgaagccaaagccaggtgtaaaagaatggacaacacagcggagagatacttggagatgaagtcaactttttcacttgatgaagaatttgattatgacaatgtgatgctaa
cctccaagttagtctgcagagatagagaacatcaagagctatgcaagcagcagaagagaaggacaccagcccagacttagagaatcctgtgactga