



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

IGF2BP3-insulin-like growth factor 2 mRNA binding protein 3 Gene

Reference: PTXBC051296

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT98; IMP-3; IMP3; KOC; KOC1; VICKZ3; insulin-like growth factor 2 mRNA-binding protein 3; IGF-II mRNA-binding protein 3; IGF2 mRNA-binding protein 3; KH domain containing protein overexpressed in cancer; VICKZ family member 3; cancer/testis antigen 98; i		

For research use only.

GENE DESCRIPTION	insulin-like growth factor 2 mRNA binding protein 3
GENE ID	10643
ORIGIN SPECIES	Human
ACCESSIONS	BC051296
VECTOR	pDONR223
NCBI SYMBOL	IGF2BP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacaaactgtatatcggaaacctcagcgagaacgcgccccctggacctagaagaatgtatctcaaggagcgcaagatcccggtgtcgggacctctctgtgtgaagactggctacgcgttcgtgtgactgcccgacgagagctgggcccctcaaggccatcgaggcgcttcagggtaa
aatagaactgcacgggaacccafagaagtgaagcactcgtccccaaaaagcaaaagatcgaaacctcagatagcaaatatccccctcattacaagtgggaagaaatggggcaaccttga