

TMEM25-transmembrane protein 25 Gene

Reference: PTXBC051841

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 25; 0610039J01Rik		
For research use only.			

GENE DESCRIPTION	transmembrane protein 25
GENE ID	84866
ORIGIN SPECIES	Human
ACCESSIONS	BC051841
VECTOR	pDONR223
NCBI SYMBOL	TMEM25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtgcctccagggcccgccctccggcacaactgctgctcctgccagccctctgagctcaggttggggggagttggagccacaataatggtcagacctgggctgagcgggcacttgggagaatgaacgccacgcttcacgtccgggtggcagggggcctgg
caccocccagattggcctggtatctggatggacagctgaggaggccagcactcaagactgctgagcgtgggaggggagggcctctctggaggcaccagcacttcactgtcactgccatcgggcccagcatgagctcaactgctctctcagggacccagaagtggccgatcagcc
aacgcctctgtcatccttaattgtcaattcaagccagagattgcccgaatggcgccaaatgaccagggaagctcagggcccgccctcctgtgttgcctgtttgccctgtgtgctgccaacccgcccgaatgtcacttggatcgaccaggtggccagtgactgtcaacacctgactt
cctgggtgtggatgctgcagaactaccctgggtcaccaaccacgggtgcaactgcaactcgcagcctggcacacaactctcgggtggggccacaatgacgtgggtgtcaccagtgcgtcgttccagccccaggccctccggcaccatctctgatatcaagtactccaaca
acctaataactcaacaacgtgcgctgccacgggagacaatgtccctccgtccaactcagctcaatgacctcactccagattccagagcagtgaaaccagcagaccggcagatggctcagaacaacagccggccagagcttctggaccggagcccgccgctcctcaccagcc
gaggtttcatccgctccagtgctgggctatatctatcgagtgccagcgtgagcagtgatgagatctggctctga