

VASH2-vasohibin 2 Gene

Reference: PTXBC051856

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	vasohibin-2; testicular tissue protein Li 222; vasohibin-like protein		
For research use only.			

GENE DESCRIPTION	vasohibin 2
GENE ID	79805
ORIGIN SPECIES	Human
ACCESSIONS	BC051856
VECTOR	pDONR223
NCBI SYMBOL	VASH2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacggctcgcggccgacactcaccgtgccccaccccaaggcgccaaaggcaccgggtccggagcagccacgcggcccgctgagcctcgccaccagcggggctcagaggaggagacaaagacgggggtgctgttccagtcacaagagcggcttcccca
tcgacagccacacctggggcgcctggtgatgcacgtggccaaggtgcaccctaagggggggagaatggtggcgccatcagggaacgcgcctcttggcgaagcctcaatccccaggtcccaaacacagggctgcatgacgatccagactggctccaggcgcagcaatta
catgaagaccctacactacttaaccaatgggcagccttcattgagcgggttcccatcagctttaaactacttctcaggaaactacttccaccagttgtgctggggatttactgcaatggccgctatggctcattgggcagcgcagggctgagctgatggacaagccattgactttcg
gactctgagtgacctcatctttgactttgaggaactcttaagaatactgacacagtcagaaggtcaagattgggctgtacgtccccatgagcctcatagcttccagccattgagtggaagcagctgtctcaacgtctcaagatgctgagggctgacataaggaggagctgga
gaaatagccaggacatgagaatgaagatcctgaaacctgcaagtgcccactctccgaccaagtgaagagccgggggaaatcctgtccccagaaggagacagcgaagcccccgaggagcctggccggcgagagaagtcgctgactgctgaaaagaaggtgctgat
ctgagcacttgaatgaagtgggtatcaaatccgaatttag