

RFPL2-ret finger protein-like 2 Gene

Reference: PTXBC051910

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RNF79; ret finger protein-like 2; RING finger protein 79; ret finger protein like 2		
For research use only.			

GENE DESCRIPTION	ret finger protein-like 2
GENE ID	10739
ORIGIN SPECIES	Human
ACCESSIONS	BC051910
VECTOR	pDONR223
NCBI SYMBOL	RFPL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgcaacttccaagaagcaagcagctgtcccgctgtcagactatctggaaaaacaaatgtccctggagtgatgcgcgctgtcctcaagtgcattaatcactgcagaaggagcccatggggaggatctactttgtgtctctccatgggtctcgggaagaacaaatca
ggcgcaatcgccagctagagagggctgtcccatcaaggaaactggagcccaagctgaagaagattctgcagatgaacccaaggatgcgggaagtccaagtgcataatgacattggatgcccaacagccaacaacttctctctattctgacgacctcaggagcgtccgaagtggg
cgcatcagacagaatcgccgaagaccttgcgagagattgacgtgtccgtttgcatctgtgggtccctcgttttacctgtgcccacactgctggaggtggacgtgggaacaagcacagaatgggacctgggagctctgcagagaatctgttcaccgcaaggaggatccagctgac
cacagagctggattctggactgtgagtttgagggatggagccgcctctctgccaccacgtgcccgtgacttctctctgtagaccgcaagttacagcagtggtgggattttctgatatgggcatgcagaacgttctcttttgatgctgaaagtgggtccatgtctatcattcaggagc
gtatctgtgaggagccattgcgccattttgtctcttcagttccacctaattggtgatcaagggtgtctgtgagcatctgtcctttgatgaactcaggcactactgatgctccagtcctgtgggagggccaataa