



Reference: PTXBC052269

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BTBD4; BZEL; RINZF; ZNF340; dJ583P15.7; dJ583P15.8; zinc finger and BTB domain-containing protein 46; BTB (POZ) domain containing 4; BTB-ZF protein expressed in effector lymphocytes; BTB/POZ domain-containing protein 4; zinc finger protein 340; zinc finge		

For research use only.

GENE DESCRIPTION	zinc finger and BTB domain containing 46
GENE ID	140685
ORIGIN SPECIES	Human
ACCESSIONS	BC052269
VECTOR	pDONR223
NCBI SYMBOL	ZBTB46
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr