

NAT11-N-acetyltransferase 11 (GCN5-related, putative) Gene

Reference: PTXBC052298

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NAT11; PATT1; N-alpha-acetyltransferase 40; N(alpha)-acetyltransferase 40, NatD catalytic subunit, homolog; N-acetyltransferase 11 (GCN5-related, putative); N-alpha acetyl transferase 40; natD catalytic subunit; N(alpha)-acetyltransferase 40, NatD catalyt		

For research use only.

GENE DESCRIPTION	N-acetyltransferase 11 (GCN5-related, putative)
GENE ID	79829
ORIGIN SPECIES	Human
ACCESSIONS	BC052298
VECTOR	pDONR223
NCBI SYMBOL	NAT11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggagaaagtcagcaaaagccaaggagaagaagcagaagcgggtggaggagcgcagcagccatggatgccgtttgtgccaagtgagacgtgccaaacaggctggagaccctctggaggctttccagtggtcaagaatatgatagaacgggtgaatgtctccattgaatgtaa
gcgagtgctggactggagccagccaccgtggatggccttgacctgacccaaacgaatatgcaaacatgatgagcagagcagtggtggctggaaagaccgagagaaacgggaggaaatgacagatgaccgagcctgtacatcgcgtgggaaaacagctccgtccctg
ttgcctttctcactccgtttgacgtggagtggtggatgaagtcctgtactgtatgaagtgcagttggaaagcaagtgccgagaaagccctggggaagttctcatagatcctgcagctcatggccaacagcacacagatgaagaaggttatgtaacagtatttaacacaatca
tggctcctaccagttctcagagaagcgttgcaattgaaattgatgactctccccagcatgtccggttctgtggggaggattgtctctatgagatcctgagccggaggaccaagttggggacagccatcctccacgcgggtgggcactgtgtgtgctgtgcactga