

FAM168A-family with sequence similarity 168, member A Gene

Reference: PTXBC052341

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM168A; KIAA0280; TCRP1; tongue cancer chemotherapy resistance-associated protein 1; family with sequence similarity 168 member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 168, member A
GENE ID	23201
ORIGIN SPECIES	Human
ACCESSIONS	BC052341
VECTOR	pDONR223
NCBI SYMBOL	FAM168A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacctgtttacagcccgctgcagcctgggctccttatggcaaccctaagaacatggcctacacgggttacccacagcctatccagcagctgcccctgcctacaatccagcctgtacccaccaatagtcaggttatgctccagcaactctgctgatgaaacaggcctggccaca
gaactcgtcttctgtggcactgaaggcaccctccacccagtgagaccgggacccgagaaacgaactacaaagcctctctgctggctttcagatatactgcgggacaccatacaaggctccaccgaccagagtaaacactgtccacccccctactccccatcaccacccctat
cagacggccatgtatccaatcagaagtgcctacccacagagaatctgtatgccaggggcctactacacacacggcgtgtatgctgccagcctcatgtcatccacataccacggctgtccagcccaacagcattccctctgtatctacccagcactgttgcgccccgaggacca
acgggtgtggccatgggcatgggtggcaggcaccaccatggcaatgtcagcaggtaccctgtgactacacccagcacacggcgattggggcacaccctgtctccatgccacatataggcccaaggaaacccctgcgtacagctactgtccccacactgttaa