

MRPL10-mitochondrial ribosomal protein L10 Gene

Reference: PTXBC052601

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L10MT; MRP-L10; MRP-L8; MRPL8; RPML8; 39S ribosomal protein L10, mitochondrial; 39S ribosomal protein L8, mitochondrial; L8mt; mitochondrial ribosomal protein L10		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L10
GENE ID	124995
ORIGIN SPECIES	Human
ACCESSIONS	BC052601
VECTOR	pDONR223
NCBI SYMBOL	MRPL10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcggcgtggcgggatgctgcgaggggtctctgccccaggcggcgctgctaccctccagactgtccgctatggctccaaggctgttaccgccaccgtctgtgatgcacttcagcggcagaagctgatggctgtgactgaatatccccccgaaccagccat
ccaccatcatgctgcatctctccagccccccacaggaggagataggctcatcaggcttctccggcgagatagcagcagtttccaggacaaccgaatgatagccgtctgcagaatgtggctctgagtcgagaggacaagcttcttatgcgacaccagctcggaacacaa
gatcctgatgaaggtcttccccaccaggtcctgaagccttctggaggattccaaagtacaaaatctgtgccccctttttgtggggcacaacatgctgctggtcagtgaaagagcccaaggtaaggagatggtaggacttaaggactgtgcatctcctgcccgtgctaggtgctgatt
gatgacaccatctcagcagggcgtttatcaactactccaagctccccagcctgccccgtgtgcagggggagctgtagaggccctcacctgcctcacagcccagaccactccctgctccagcaccagccccctccagctgaccacctgttgaccagtacatcagagagcaacg
cgagaaggattctgcatgtcgccaatgggaagccagatcctgacactgttccggactcgtag