

LOC728131-ribosomal protein S6 pseudogene Gene

Reference: PTXBC052613

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	ribosomal protein S6 pseudogene
GENE ID	728131
ORIGIN SPECIES	Human
ACCESSIONS	BC052613
VECTOR	pDONR223
NCBI SYMBOL	LOC728131
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagctgaacagctccttctagccattggctgccagaaactcattaaagtgaagatgaacgcaaacttcgtactttttatgagaagcatatggccacagaaattgctgttatgctctgggtgaagaatggaagggttatgtgtccgaatcagtggtgggaacaacaagaaggtttctt
gtttgttggtgggaacaacaagcagggtatcttgacccatggccatgccgcctgctgttacattcctgttatatatagaccaaggaaagctggagaagaagagaaaattag