



Reference: PTXBC053348

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OAT4L; RST; URAT1; solute carrier family 22 member 12; organic anion transporter 4-like protein; renal-specific transporter; solute carrier family 22 (organic anion/cation transporter), member 12; solute carrier family 22 (organic anion/urate transporter)		

For research use only.

GENE DESCRIPTION	solute carrier family 22 (organic anion/urate transporter), member 12
GENE ID	116085
ORIGIN SPECIES	Human
ACCESSIONS	BC053348
VECTOR	pDONR223
NCBI SYMBOL	SLC22A12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr