



Reference: PTXBC053990

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAPP2; pleckstrin homology domain-containing family A member 8; PH domain-containing family A member 8; phosphatidylinositol-four-phosphate adapter protein 2; phosphoinositol 4-phosphate adapter protein 2; pleckstrin homology domain containing, family A (		

For research use only.

GENE DESCRIPTION	pleckstrin homology domain containing, family A (phosphoinositide binding specific) member 8
GENE ID	84725
ORIGIN SPECIES	Human
ACCESSIONS	BC053990
VECTOR	pDONR223
NCBI SYMBOL	PLEKHA8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1