

SPANXC-SPANX family, member C Gene

Reference: PTXBC054023

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT11.3; SPANX-C; sperm protein associated with the nucleus on the X chromosome C; SPAN-Xc protein; cancer-testis-associated protein CTp11; cancer/testis antigen 11.3; cancer/testis antigen family 11, member 3; cancer/testis-associated protein of 11 kD; nu		

For research use only.

GENE DESCRIPTION	SPANX family, member C
GENE ID	64663
ORIGIN SPECIES	Human
ACCESSIONS	BC054023
VECTOR	pDONR223
NCBI SYMBOL	SPANXC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacaacaatccagtgcggcggtgaagaggagctccctgtgaatccaacgaggtgaatgagacgatgccggagacccaactgggactcagaccgcaacctgctctaaaaaatgaaacatctgagtcctcgaccatactagtgttcgctacaggaggaacgtg
aaaagaacatctccagagaactgctgaatgaccacgcccagagagaacagaatcaacccctccaaatggaggaggaggaattcatggaataatggttgaaatacctgcaagtag