

A1CF-APOBEC1 complementation factor Gene

Reference: PTXBC054873

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACF; ACF64; ACF65; APOBEC1CF; ASP; APOBEC1 complementation factor; APOBEC-1 stimulating protein; apo-B RNA editing protein; apobec-1 complementation factor (ACF) (ASP)		

For research use only.

GENE DESCRIPTION	APOBEC1 complementation factor
GENE ID	29974
ORIGIN SPECIES	Human
ACCESSIONS	BC054873
VECTOR	pDONR223
NCBI SYMBOL	A1CF
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaatcaaatcacaatcggggatggattgagcggcactcagaagggaagcagccctccgcgactgtccagcgcacaggatagcttggccaggaaaatggacaagaaaatgatggtggccctccacctggtgggatgctgcacccctgaaaggggctgtgaaattttatt
ggaaaacttccccgagaccttttggagtgagcttataccattatgtgaaaaaatcgtaaaatttatgaaatgagaatgatgatggttttaafggcaacaatagaggatgcatgttaacattttcaaataaagtggagccaagaatgcaatcaagcaacttaataattatgaaattaggttaa