

TINAG-tubulointerstitial nephritis antigen Gene

Reference: PTXBC056235

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TIN-AG; tubulointerstitial nephritis antigen		
For research use only.			

GENE DESCRIPTION	tubulointerstitial nephritis antigen
GENE ID	27283
ORIGIN SPECIES	Human
ACCESSIONS	BC056235
VECTOR	pDONR223
NCBI SYMBOL	TINAG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggaccggatataagatcttaactcttcttactacagaaatctggatggagaagcagttatctcaagagaaagtgacctagaggcttattcactaggaaacacaccgtttgcaagggtactcgattcaaaagagccatttccaaggcgaactgtagaattttggctgttg
aagacagagatgatggctgtgctactgagttctatgcggcgaatgcgtgtgctactgtgataaattctgtgacagagaaaattctgattgctgctgactacaagtcctttgcccgtgaagagaaagaaatggcctcctcacacacagccttgatccagaagggtgcttcaagatggcacaac
attatgaagaggatcagtaattaaagaaaactgcaactcctggtaa