

LSM14B-LSM14B, SCD6 homolog B (*S. cerevisiae*) Gene

Reference: PTXBC057387

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LSM14B, SCD6 homolog B; C20orf40; FAM61B; FT005; LSM13; RAP55B; bA11M20.3; protein LSM14 homolog B; LSM14 homolog B; RNA-associated protein 55B; family with sequence similarity 61, member B; hRAP55B; LSM family member 14B		

For research use only.

GENE DESCRIPTION	LSM14B, SCD6 homolog B (<i>S. cerevisiae</i>)
GENE ID	149986
ORIGIN SPECIES	Human
ACCESSIONS	BC057387
VECTOR	pDONR223
NCBI SYMBOL	LSM14B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcggctcctcagcaccctgtatctgggcagcaagatcagctcatctccaaggcgcagatccgctacgagggcattctctacaccatcgacaccgacaactccaccgtggcgctcgccaaagtgggtcctttggcactgaagaccgtcccacagataggcctgcgccccag
agaggagatttatgagtacatcatcttcgaggaagtgacatcaaggatcactgtgtgtgaacctccgaaagctcagcacactcccgaggatcccgccattgttcagttctccctgggttctgctccgctcgcccttcagccgcagctgcttacagcccttccgaggatggcg
ccctacggcccgctggcgccagctccctgctcagccagcagatgcccctccctgggtctagagaagttgtaagccctccagcctcagccgctgttccagccctagctcttctccatctccacagcccggttcagagcttgactgtctcagagccacagcagctcactccaaggg
taacagcagtttgggagaactgcagctgtcttgcacacattctgagagcaggggtaagcagcagacagaatgacagtagctgtagcagatcatttccaagccatgcagccgtga