

GJB7-gap junction protein, beta 7, 25kDa Gene

Reference: PTXBC057813

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CX25; bA136M9.1; connexin25; gap junction beta-7 protein; connexin-25; gap junction protein, beta 7, 25kDa; gap junction protein beta 7		

For research use only.

GENE DESCRIPTION	gap junction protein, beta 7, 25kDa
GENE ID	375519
ORIGIN SPECIES	Human
ACCESSIONS	BC057813
VECTOR	pDONR223
NCBI SYMBOL	GJB7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttggatgttcctcagagatcctcctgagtgagtaataataactccactgggactggatgttgctggctgtcgtgtttgtctccgtttgctgtctacatgggtggcagcagagcagctgttggaagatgagcagaagagttgagtgcaacagtagacagccgggtgcaaaa
atgtgtgtttgatgactcttccccattcccaagtcagactttggccttacaactgataatggctccacacctcactctgtgtgttttacatgtagcctatcatgagggtagagagaaaaggcacagaagaaactctatgtcagccagggtacaatggatggggcctatgtacgcttat
cttatcagcctcattgttaaaactggtttgaaattggctcctgttttattttataagctatatgatggcctttagtgttccctacctataaagtgtgattgaagcctgtgcccaacactgtggactgtctcatcctcaaaaccactgagaagacgatcttcatctcttctgtgtatcacctcatgtgtg
tattgtgtgaatttcattgaactgagttttgttctcaagtgtttattaagtgtgtctccaaaaatatttaaaaaaacctcaagtcctcagtggtga