

LRRC57-leucine rich repeat containing 57 Gene

Reference: PTXBC058935

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	leucine-rich repeat-containing protein 57; leucine rich repeat containing 57		
For research use only.			

GENE DESCRIPTION	leucine rich repeat containing 57
GENE ID	255252
ORIGIN SPECIES	Human
ACCESSIONS	BC058935
VECTOR	pDONR223
NCBI SYMBOL	LRRC57
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaacagtcgcctccgctcatgtggaacggcgagaaactggtgtctttcagcttaaggaccgagggctgaccgagttcccgagactgcagaagctgacgagcaatctcaggaccatcgactgtccaaacaagaatcgaaagcctaccgctttgctgataggaaa
gttcactctgtgaagagcctctccctgaacaacaactgactgttctgctgatgagatgcaatctgaaaaactagagacgctaagcctaacaacaatcacttagagagctgcccgtctacctttggcgaactctctgcccctcaagacctgagcctctctgggaaccaactggga
gcattacctccccactttgtagcctacggcactggatgfatggaatctcttaagaaccagattcgaagtatactgactcagtgaggagagctgcaagtcacgaactcaactcaaccagaatcagatactcagatctcagtgagatattctgtccacgecttaaaattctgcctgg
aagagaattgtctgagctcagcatgttccccagagcactcagtgattccccagatctgtctgctgtggaaggcaatcttttgaataaagaaactcgagaactggaaggctatgataagtacatggagaggttcacagccaccaagaagaagtttgcgtga