



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CREB5-cAMP responsive element binding protein 5 Gene

Reference: PTXBC059400

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRE-BPA; CREB-5; cyclic AMP-responsive element-binding protein 5; cAMP response element-binding protein CRE-BPA; cAMP responsive element binding protein 5		

For research use only.

GENE DESCRIPTION	cAMP responsive element binding protein 5
GENE ID	9586
ORIGIN SPECIES	Human
ACCESSIONS	BC059400
VECTOR	pDONR223
NCBI SYMBOL	CREB5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccagctccatgctgggacccctgcccaacctcaaatgcaggatgtctccgcctgttgtagcgaattggagcgacaataatgcagtgaaactccagatcatatgggagatcaaggctccaaatctcagcaaccctctgctgtctcccaggttcacccaatgcattcagaagccaaaatga
cgttgaaagctgcatgtgactacacccctgtcgcctatgcaattgggaacatgaacacatggagacatcatgagatgaggtggctccgggacagcagacgccacatcatatgcaatgcacccgatcagcaccagcagctacccacatcccttaccacacagca
ccggccacccagcagcattcattcccttccctacccacagcaggaacatcattaccactccattccacatctatgcacacccagatcaccacagctcccatctgcctctcagcgcacgaacccagcagctgtttccacgaacaacaaagatgcagctccgcaac
cagacaatacagccacccccagccacagggggggccggcggaaggggtgtgtagacgagatccggacgagaggccggcgaaattcttgaacggaaacggggcagctgccaccctgtgcagacagaagaggaaaggtctgtgtgatgtcattggaaaagaagcagaagaactca
cccgacaacaatctgcagcttcagaataagttgtctctgttgaataatgagctggccagctcaaacaggtctgtttgaacataaagaatgcgccataacagccgatcgaagaagaatcacaaggatattcaagctccagagagtagccctctctgactgtctgcccaactgtctccagcaac
aaatcagctacataatcagactacttctctctctgctgacgcagagctccacctcagccacgtccaccatcacaagaacacactgaatccgattcttaa

More info: Please contact us for any question: contact@proteogenix.fr