

LOC649946-ribosomal protein L23a pseudogene Gene

Reference: PTXBC060042

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	ribosomal protein L23a pseudogene
GENE ID	649946
ORIGIN SPECIES	Human
ACCESSIONS	BC060042
VECTOR	pDONR223
NCBI SYMBOL	LOC649946
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagaagacagaagacaacacgcactgtgtcactgtgaatgttaaagccaacaagcaccagaagaagctctatgacactgatgtggccaaggtcaacacctgattcggcctgatggagagaaggaacatgttcgactggctcctgattacaatgctttggatgttgccaacaaat
agggatcatctcaactgagtcagttggctaattctaaatatatgtatatctttcaccattaa