

HIST2H2AC-histone cluster 2, H2ac Gene

Reference: PTXBC060324

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2A; H2A-GL101; H2A/q; H2AFQ; histone H2A type 2-C; H2A histone family, member Q; histone 2, H2ac; histone H2A-GL101; histone H2A/q; histone IIa; histone cluster 2, H2ac; histone cluster 2 H2A family member c		

For research use only.

GENE DESCRIPTION	histone cluster 2, H2ac
GENE ID	8338
ORIGIN SPECIES	Human
ACCESSIONS	BC060324
VECTOR	pDONR223
NCBI SYMBOL	HIST2H2AC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgtgtgtgcaacaaggaggcaagcccgcgccaaggccaagtcgctcgtcccgctggtccctcagttcccgtagggcgagtgaccgcttgcgcgcaaaaggcaactacgaggagcgggtggggccggcgcccggtctacatggcggtctcgtgagtagc
tgaccgagatcctgagctggcgggcaacgcgctcgggacaacaagaagacgcgcacatccctcgtcactccagctggccatccgcaacgacgaggaactgaacaagctgctgggcaaaagtcaccatgccaggcggtttttgcctaaccagggccgttctgttac
caagaaaaccgaaagccacaagccaaaagcaataa