



Reference: PTXBC060818

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AP180; CALM; clathrin coat assembly protein AP180; 91 kDa synaptosomal-associated protein; assembly protein, 180kDa; clathrin coat-associated protein AP180; phosphoprotein F1-20; synaptosomal-associated protein, 91kDa homolog; synaptosome associated prote		

For research use only.

GENE DESCRIPTION	synaptosomal-associated protein, 91kDa homolog (mouse)
GENE ID	9892
ORIGIN SPECIES	Human
ACCESSIONS	BC060818
VECTOR	pDONR223
NCBI SYMBOL	SNAP91
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]