

RPL13A-ribosomal protein L13a Gene

Reference: PTXBC062537

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L13A; TSTA1; 60S ribosomal protein L13a; 23 kDa highly basic protein; tissue specific transplantation antigen 1; ribosomal protein L13a		

For research use only.

GENE DESCRIPTION	ribosomal protein L13a
GENE ID	23521
ORIGIN SPECIES	Human
ACCESSIONS	BC062537
VECTOR	pDONR223
NCBI SYMBOL	RPL13A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggaggtgcaggtcctggtgctgatggtcaggccatctcctggccgcctggcgccatcgtggttaaacaggtactctggggcgggaaggtggtgctgtacgtgtgaaggcatcaattcttgcaatttctacagaacaagttgaagtacgtggttctcctcgaagc
ggatgaacaccaacccctcccgagggccctaccacttccgggccccagccgcattcttggcgaccgtgcgaggtatgctgccccacaaaacaaagcgagggccgctctggaccgtctcaaggtgtttgacggcatccaccgccctacgacaagaaaaagcggatggtg
ttcctgctgccctcaaggctgtgctttaaagcctacaagaaatttgcctatctggggcgcctgggtcagcaggttggctggaagtaccaggcagtgacagccaccctggaggagaagaggaaagagaagcaagatccactaccggaagaagaacagctcatgaggctacggaa
acaggccgagaagaacgtggagaagaaattgacaatacacagaggtcctcaagaccacggactcctggtctga