

## CDH26-cadherin-like 26 Gene

Reference: PTXBC062570

|                |                                                                                           |                     |                |
|----------------|-------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                       | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                      | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | VR20; cadherin-like protein 26; cadherin-like 26; cadherin-like protein VR20; cadherin 26 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | cadherin-like 26                                                                                  |
| <b>GENE ID</b>          | 60437                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC062570                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | CDH26                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgaaacctgtgatggacatgttcagatgttgaaggccagaggccggctctgtctatctgcacagctgcagcaggaccacgcaggagtgtaaggctaccagatgccacaatgcacagacaactcctggctccggtggaaggagggatggcagagacattgaatcagaaactcca  
tggtgccaatgtgctggaagatgaccccggtacctacctcacgtctacagcagggaaggaggagtggtgaggggcccatccctcagctctctgcccagcttggaacaggagttgcaacctgattgctggaactctttgggttcaaaagcgactccgttgaggaaatattcagagtcagg  
tgttcctcctaa