

PRO0461-hypothetical LOC652276 Gene

Reference: PTXBC062576

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRO0461; PDK1; PDPK2; PDPK2P; 3-phosphoinositide-dependent protein kinase 1; 3-phosphoinositide-dependent protein kinase 2 pseudogene; PkB kinase like gene 1; 3-phosphoinositide dependent protein kinase 1		

For research use only.

GENE DESCRIPTION	hypothetical LOC652276
GENE ID	652276
ORIGIN SPECIES	Human
ACCESSIONS	BC062576
VECTOR	pDONR223
NCBI SYMBOL	PRO0461
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggggaatcactgcgagctcctgccgtggccgtggcaggtcggggcgggttgggtggcttcttggcctccctaaagcgcggggtcagcgtcctggccagcggccagcagcaggtccaagtgggtccggtctacagcggcgccactactctcaccactgg
gcagacgcgcgtgtgtcaggaccgaaatcctctgtacctcttgagccagccgaccccaaccggactcggacaagacggagtgtttgtctgttggccaagctggagtacaatggcacaatcttggtcaccacaacctctgccactgggttcaagcaggtctctctcagtcctct
ga