

PFDN5-prefoldin subunit 5 Gene

Reference: PTXBC062671

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MM-1; MM1; PFD5; prefoldin subunit 5; c-myc binding protein; myc modulator-1		

For research use only.

GENE DESCRIPTION	prefoldin subunit 5
GENE ID	5204
ORIGIN SPECIES	Human
ACCESSIONS	BC062671
VECTOR	pDONR223
NCBI SYMBOL	PFDN5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcagctctattaacatcacggagctgaatctgccgcagctagaatgctcaagaaccagctggaccaggaagtggagttctgtccacgtccattgctcagctcaaaagtgtacagaccaagtatgtggaagccaaggactgtctgaacgtgctgaacaagagcaacgaggggaaa
gaattactcgtccactgacgagttctatgtatgtccctgggaagctgcatgatgtggaacacgtgctcatgcatgtgggaactgggtactatgtagagaagacagctgaggatccaaggacttctcaagaggaagatagatttctaaccaagcagatggagaaaatccaaccagctcttc
aggagaagcagcgcctatgaacagccgctcatggaatgatgagtcagaagattcagcagctcacagccctgggggcagctcagggtactgctaaggcctga