

CST2-cystatin SA Gene

Reference: PTXBC062679

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	cystatin-SA; cystatin 2; cystatin S5; cysteine-proteinase inhibitor; salivary cysteine (thiol) protease inhibitor; cystatin SA		

For research use only.

GENE DESCRIPTION	cystatin SA
GENE ID	1470
ORIGIN SPECIES	Human
ACCESSIONS	BC062679
VECTOR	pDONR223
NCBI SYMBOL	CST2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctggccctgtgcacctgctcctgtgcccacccagcgtgtggccctggcctggagccccaggaggaggacaggataatcagggtggcatctatgatgcagacctcaatgatgagcgggtacagcgtgccctcactttgtcatcagcgagtataacaaggccactgaa
gatgagtactacagacgcctgtcgggtgtctacagccaggaggagcagatcgtggcgggggtgaattacttcttgacatagagggtggccgaacctatgtaccaagtcagcccaactggacacctgtgccttccatgaacagccagaactgcagaagaacagttgtctcttcc
agatctacgaagttccctgggaggacagaatgtccctgtgaattccaggtgtcaagaagcctag