

XAGE3-X antigen family, member 3 Gene

Reference: PTXBC062680

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT12.3a; CT12.3b; GAGED4; PLAC6; XAGE-3; pp9012; X antigen family member 3; CT12.3; G antigen, family D, 4; cancer/testis antigen 12.3; cancer/testis antigen family 12, member 3a; cancer/testis antigen family 12, member 3b; g antigen family D member 4; pl		

For research use only.

GENE DESCRIPTION	X antigen family, member 3
GENE ID	170626
ORIGIN SPECIES	Human
ACCESSIONS	BC062680
VECTOR	pDONR223
NCBI SYMBOL	XAGE3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgattggcgaggagatcaacatagccctagccgaggagagaagtgtaccacctctgagctgattggcctatgctggagcccggtgatgaggagcctcagcaagagggaaccaccaactgaaagtcgggacacctggtcaggagagagaagaagatcagggtgcagct
gagactcaagtgcctgacctggaagctgactccaggagctgtctcagtcagcaagactgggggtgaatgtggaatggtcctgatgaccaggggaagattctgccccaaatcagaacaatttaaatgccagaaggaggtgacagggaaccacaggtttaa