

C6orf58-chromosome 6 open reading frame 58 Gene

Reference: PTXBC062712

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0762 protein C6orf58; LEG1; protein LEG1 homolog; liver enriched gene 1 homolog; chromosome 6 open reading frame 58		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 58
GENE ID	352999
ORIGIN SPECIES	Human
ACCESSIONS	BC062712
VECTOR	pDONR223
NCBI SYMBOL	C6orf58
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggatgtataaatcatattgaatcagacagccaggtatttgcacaaattgcaccagataatgaacagaattttatgggggttgcctctgcagtagtgctggcaatagtagcagagcagattagctgatccaacccgaaggacaaactgtggctatgaatctggagatcatatgtgcac
tctgtggacagttgtgggctgattgaattatttctgtcttcattacccttcttgcgtcggttgattctggtgtaattgggagatcatcagaccaagtcaggcttttgcacccacccaagaatgagagggaagttttgtatgatgtttctagctgtcttcctccctgagacaatgaacaagtgg
aacacctttaccagatttgcagtcaccttttagtaagtttgatgatctgttgagtaactatgggctgcacacacttcaaccttggcagataatacaaaagttttgaagacagatagattatttctaaagcagaagcgcatttgagagaagttgggtactggctgtggatcatttagctgcagtc
ctcttctacaaccttgattagatcatataagttccagaaggcgatgccaccagaaattcttctaactgatgtagcccttcatcagtgactttactgctttcagaatgtagtcctggttctctaaatatgcttgacaatgtggataaatctataggttatctttgtacagaaaaatctaattgtatata
gagatcattcggaatctagctctagaagttatggaataaactcctga