

HEATR5A-HEAT repeat containing 5A Gene

Reference: PTXBC062720

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf125; HEAT repeat-containing protein 5A; HEAT repeat containing 5A		

For research use only.

GENE DESCRIPTION	HEAT repeat containing 5A
GENE ID	25938
ORIGIN SPECIES	Human
ACCESSIONS	BC062720
VECTOR	pDONR223
NCBI SYMBOL	HEATR5A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaatagctcatagcttactgctgaatgaagaagcacaatcaactaggtgaagtcagaaggcagagttatttttgagtggttgagatactggagaagctcttggcaaccagcaggaatgatgaaggagagaacagaagactctgttgaacagctcctgtctttgtgaacagct
ccccaggcctcctaccgcaaacctgtgctaagaatctagccatactttatagtattggagacacattctccgttcaggaagcaatcgataaatgtaatgatcttattcgtagcaaatgatctccaagtattcttcccactaagctgtaa