

**BOLA2-bolA homolog 2 (E. coli) Gene**

Reference: PTXBC062756

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BOLA2A; BOLA2B; My016; bolA-like protein 2; BolA-like protein 2 member A; bolA homolog 2; bolA family member 2		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	bolA homolog 2 (E. coli)
<b>GENE ID</b>	552900
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC062756
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	BOLA2
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

atggaaactcagcggccgaatactccgcgagaagctgcagcgggacctggagcggagcatgtggaggtggaggacacgacctcaaccgttgctcctgtagctccgagtcctggtggtgtcggccaagttcgagggaaccgctgcttcagagacacaggttctgtacagagtaa