

FAM13A1-family with sequence similarity 13, member A1 Gene

Reference: PTXBC063126

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAM13A1; ARHGAP48; protein FAM13A; FAM13A1_v2 protein; family with sequence similarity 13, member A1; family with sequence similarity 13 member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 13, member A1
GENE ID	10144
ORIGIN SPECIES	Human
ACCESSIONS	BC063126
VECTOR	pDONR223
NCBI SYMBOL	FAM13A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggcaggagctctagccatctgtcaaaagcagcggctggctgaaagaagacatgaaaaagatagtgccagtgccattaaatgaacagaaggattttacattcagaagttatgtggagtcagtcctcaagaactgaacggcaggggctcaccgagaatggcattccagcag
tagtgggaatagtggaatatttgacgcagcatggacttaccgaagaaggctcttttaggggtgaatggaacgtgaagggtgggaacaacttcgactgaagttcgagagtggaagtgcccgaggagctcgggaaggacgggtgatgtctgctcagcagccagctctgttgaagctgtttctga
gggagctgccctgacagctgatcacctcagcgttgcagcctcattcattcaactcttcaggatggcagaatgatgttcaggagagtagcttaagagacttaataaaagagctgccagacaccactactgcctcctcaagtactttgccagttcttgacaaaagtagccaagcatcatgtg
cagaatcgcataatgttcacaatctgccactgtatttgggccaattgtttcagtaa