

ZNF468-zinc finger protein 468 Gene

Reference: PTXBC063143

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein ZNF468; zinc finger protein 468		

For research use only.

GENE DESCRIPTION	zinc finger protein 468
GENE ID	90333
ORIGIN SPECIES	Human
ACCESSIONS	BC063143
VECTOR	pDONR223
NCBI SYMBOL	ZNF468
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatggcgaaggttttaatacaaaaagcaaccttgcacatcaccatagactcatactggagagaaacctacaatgggagggaatfgacaaaagtttcagttgtaaatcaaaccttaaacacataagaaaatcacattgaagagaaccatacagaggttaaggtttgacaaggttttcataataatgcataatctggcaaacatactagaattcacactggagagaaactcattataagttaatgagtggtggcaagaccttggcaaaattcacacctgtaa