



Reference: PTXBC063283

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ADAMTS21; KNO2; MMCAT; A disintegrin and metalloproteinase with thrombospondin motifs 18; a disintegrin-like and metalloprotease (repolysin type) with thrombospondin type 1 motif, 18; a disintegrin-like and metalloprotease (repolysin type) with thrombos		

For research use only.

GENE DESCRIPTION	ADAM metallopeptidase with thrombospondin type 1 motif, 18
GENE ID	170692
ORIGIN SPECIES	Human
ACCESSIONS	BC063283
VECTOR	pDONR223
NCBI SYMBOL	ADAMTS18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]